

Banked hBE Vials by Genotype

| Genotype | #Vials | #Patient Samples |
|------------------------|--------|------------------|
| F508del/F508del | 479 | 9 |
| F508del/R553X | 415 | 4 |
| WT/WT | 293 | 3 |
| Unknown | 284 | 3 |
| F508del/1811+1.6kbA->G | 213 | 1 |
| K710X/L467P | 192 | 1 |
| F508del/3120+1G->A | 173 | 1 |
| F508del/W1282X | 166 | 2 |
| F508del/3140-26A->G | 159 | 1 |
| F508del/N1303K | 132 | 3 |
| F508del/A559T | 131 | 1 |
| F508del/E1371X | 131 | 1 |
| F508del/3849+10KbC->T | 110 | 1 |
| F508del/V520F | 110 | 1 |
| F508del/1898+1G->A | 104 | 1 |
| F508del/1154insTC | 101 | 1 |
| F508del/296+1G->T | 95 | 1 |
| R553X/W1282X | 93 | 1 |
| F508del/3272-26A->G | 89 | 1 |
| F508del/E60X | 84 | 1 |
| F508del/G542X | 65 | 2 |
| F508del/M1101K | 64 | 1 |
| F508del/G551D | 62 | 1 |
| F508del/1078delT | 57 | 1 |

| | | |
|-------------------|----|---|
| W1282X/W1282X | 55 | 1 |
| F508del/621+1G->T | 54 | 1 |
| WT/R792G | 54 | 1 |
| G542X/621+1G->T | 54 | 1 |
| R553X/G85E | 50 | 1 |
| F508del/S549N | 46 | 1 |
| F508del/2184delA | 42 | 1 |
| F508del/Q493X | 15 | 1 |
| Y122X/Y122X | 14 | 1 |
| F508del/R334W | 11 | 1 |
| F508del/R347P | 11 | 1 |